

How hearing instruments can help

Today's technologies deliver clarity and convenience

Automatic program

Seamlessly adjusts hearing instrument settings to meet your listening needs when the environment changes.

SmartFocus™

Ability to emphasize clarity of sound in front of you or comfort of sounds around you.

Self learning

Gradually learns your situational listening preferences and applies them automatically.

LearnNow™

Instantly teach your hearing instruments your preferences for loudness, comfort and clarity across many situations.

Pinna Effect

Recreates the natural ability of the ear to collect sound.

Binaural Phone

Automatically hear the phone conversation in both ears to improve clarity and understanding.

Natural Sound Balance

Continuously monitors and adjusts the hearing instrument settings to preserve natural sound quality.

DuoLink

Both hearing instruments are adjusted when program or volume control changes are made on one side.

Directional microphone*

Improves speech understanding in noise by focusing on sounds in front of you and decreasing background noise from other directions.

MyMusic™

Enhanced music listening experience based on your preferences.

Feedback manager

Instantly suppresses annoying whistling or feedback.

Speech enhancement

Boosts speech in everyday situations to improve clarity and ease of listening.

Noise reduction

Decreases detected background noise to provide listening comfort.

AntiShock™

Instantly suppresses loud, irritating sounds without affecting sound quality.

Wind noise manager

Reduces wind noise to improve listening comfort.

Data logging

Tracks real-world usage to guide fitter adjustments and provide more effective follow-up appointments.

IntelliVent

Customized venting process to improve comfort and enhance sound quality.

Easy-t

Automatically changes to telephone program.

Easy-DAI*

Automatically changes to direct audio input program when sound is received via a DAI attachment (e.g, FM signal).

* Availability depends on the style of instrument selected. Talk to your hearing healthcare professional for more information.

Accessories

Smart Control

Discrete, easy-to-use remote control allows adjustment of hearing instrument settings.

uDirect™

Convenient hands-free interaction with your cell phone and other audio devices.

uTV™

Sends sound wirelessly from your TV to your hearing instruments.