

Apps feature overview.

Remote support, audio streaming and remote control at the tip of your clients' fingers: Most Signia hearing aids are compatible with useful apps for Android and iOS smartphones. This overview shows you the features of each app in detail.

	Remote support	Audio streaming and remote control	Remote control	Audio streaming and remote control
	 myHearing App	 myControl App	 touchControl App	 easyTek App ^A
Availability	Apple App Store/ Google Play Store	Apple App Store/ Google Play Store	Apple App Store/ Google Play Store	Apple App Store/ Google Play Store
Use case	Signia TeleCare	Hearing aids with Bluetooth connectivity	All primax ^F , all binax, Orion 2, Sirion 2, Intuis 3 and Prompt hearing aids	easyTek required
Remote control features				
Transmission signal	Bluetooth High frequency coded signal	Bluetooth High frequency coded signal	High frequency coded signal	Bluetooth High frequency coded signal
Volume adjustment	●	●	●	●
Independent volume adjustment (R/L)	-	● ^B	-	-
Program change	●	●	●	●
Mute	●	●	●	●
SoundBalance adjustment	●	●	●	●
Adjustable system sounds	-	● ^C	-	-
Spatial Configurator adjustment	● ^D	● ^D	● ^D	● ^D
Automatic listening situations	-	● ^E	-	-
Power-on-Delay configuration	-	● ^E	-	-
Battery status read-out	-	●	● ^F	●
Low battery level	-	●	-	-
Voice activity tracking	-	● ^E	-	-
Sound exposure	● ^G	●	-	-
Wearing time	● ^G	●	-	-
User guide	●	●	-	-
Adaptive streaming volume	-	● ^{E, H}	-	-
Smartphone motion sensors usage	● ^G	●	-	-
Hearing aid airplane mode	-	●	-	-
Hearing support	● ^G	●	-	-
Signia TeleCare features				
Usage statistics	● ^G	-	-	-
Daily satisfaction rating	● ^G	-	-	-
Hearing lessons	●	-	-	-
CareChat	●	-	-	-
Private label option	● ^I	-	-	-

(A) Setting in Connexx needed.

(B) Configurable via Connexx.

(C) Volume level adjustment not available with Pure 13 BT primax.

(D) Available with performance levels 7 and 5, only.

(E) Not available with Pure 13 BT primax.

(F) Output via the number of signal tones in the hearing aid.

(G) For Pure 312 Nx, Pure 13 Nx, Motion 13 Nx and Pure 13 BT primax only.

(H) Not available in the app version on the Google Play Store.

(I) After activation via customer service.