

Signia Nx product overview.



	SLIM-RIC		RIC				BTE			CLICK ITE	Custom
	Styletto Connect	Styletto	Pure™ 10 Nx	Pure™ 312 Nx	Pure™ 13 Nx	Pure™ Charge&Go Nx	Motion™ Charge&Go Nx	Motion™ 13 Nx	Motion™ 13P Nx	Silk™ Nx	Insio™ Nx IIC/CIC / ITC/ITE
Battery size	–	–	10	312	13	–	–	13	13	10	10/312
Li-ion rechargeable	✓	✓	–	–	–	✓	✓	–	–	–	–
Tcoil option (battery door)	–	–	–	–	✓	–	–	✓	✓	–	–
Direct streaming	✓	–	–	✓	✓	✓	✓	✓	✓	–	✓*
Telecare ready	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
StreamLine TV/Mic compatible	✓	–	–	✓	✓	✓	✓	✓	✓	–	✓*
myControl App	✓	–	–	✓	✓	✓	✓	✓	✓	–	✓*
touchControl App™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
miniPocket™	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Wireless programmer	NoahLink Wireless	ConnexxAir	–	NoahLink Wireless	NoahLink Wireless	NoahLink Wireless	NoahLink Wireless	NoahLink Wireless	NoahLink Wireless	ConnexxAir	ConnexxAir /NoahLink Wireless*
HiPro	–	–	✓	✓	✓	–	–	✓	✓	✓	✓

* Insio Nx ITC/ITE only

For all general enquiries or to place an order please call **01293 423700** or email **info-hearingaids.uk@sivantos.com**

Our game-changing new platform Signia Nx delivers the most natural own voice and the best speech understanding in noise, while combining the most advanced connectivity systems.

Features	Information	Performance levels		
Audiology		7Nx	5Nx	3Nx
Signal processing and Gain/MPO	Channels/Handles	48/20	32/16	24/12
OVP™ (Own Voice Processing)	Own Voice Processing (OVP™) utilizes real-time recognition of the wearer's voice to deliver a natural own voice impression.	●●●●●	●●●●●	●●●●●
3D Classifier	Unlike conventional classification systems which rely on acoustic data, 3D Classifier analyzes multi-dimensional information from the wearer's environment and communication behavior. Information on the wearer's voice activity, motion state, and acoustic data are used to tailor the hearing aid's response to the wearer's needs.	●●●●●	●●●●●	●●●
Sound Clarity	Sound Clarity provides lifelike sound which covers the full spectrum of wearing situations. Exceptional clarity is achieved with features such as: HD Spatial, extended dynamic range, extended bandwidth, EchoShield, HD Music, eWindScreen™, speech and noise management and feedback cancellation.	●●●●●	●●●●●	●●●
Speech Quality	Speech Quality delivers improved understanding of other people's voices in noise by using an exclusive technology called Ultra HD e2e. Features which utilize Ultra HD e2e include narrow directionality, Spatial SpeechFocus™ and TwinPhone.	●●●●●	●●●●●	●●●
Connectivity				
Direct Streaming	Signia Nx hearing aids can be connected directly to Apple devices for phone calls and audio streaming - without additional accessories. With the StreamLine Mic accessory, direct streaming is possible also from Android devices. With the StreamLine TV accessory, high-quality stereo sound can be streamed directly to Nx hearing aids.	●●●●●	●●●●●	●●●●●
You Adjust	Technically inclined wearers often appreciate having extra control over their hearing aids. Signia Nx hearing aids offer these wearers a comprehensive range of adjustment options via the myControl App.	●●●●●	●●●●●	●●●●
Remote hearing care				
TeleCare	TeleCare provides a structured communication approach for hearing care professionals and their clients, which is especially important during the initial stages of a fitting. In addition, TeleCare provides full live remote tuning to assist with rapid and effective troubleshooting.	●●●●●	●●●●●	●●●●●

●●●●● = highest feature performance