

Oticon miniRITE R: Oticon More™








	Oticon More 1	Oticon More 2	Oticon More 3
Speech Understanding			
MoreSound Intelligence™	Level 1	Level 2	Level 3
Environment configuration	5 Options	5 Options	3 Options
Virtual Outer Ear	3 Configurations	1 Configuration	1 Configuration
Spatial Balancer	100%	60%	60%
Neural Noise Suppression, Difficult / Easy	10 dB / 4 dB	6 dB / 2 dB	6 dB / 0 dB
Sound Enhancer	3 Configurations	2 Configurations	1 Configuration
MoreSound Amplifier™	•	•	•
Feedback Prevention	MoreSound Optimizer™ & Feedback shield	MoreSound Optimizer™ & Feedback shield	MoreSound Optimizer™ & Feedback shield
Spatial Sound™	4 Estimators	2 Estimators	2 Estimators
Soft Speech Booster	•	•	•
Frequency lowering	Speech Rescue™	Speech Rescue™	Speech Rescue™
Sound Quality			
Clear Dynamics	•	•	-
Better-Ear Priority	•	•	-
Fitting Bandwidth	10 kHz	8 kHz	8 kHz
Bass Boost (streaming)	•	•	•
Processing Channels	64	48	48
Listening Comfort			
Transient Noise Management	4 configurations	3 configurations	3 configurations
Wind Noise Management	•	•	•
Personalisation & Optimising Fitting			
Fitting Bands	24	20	18
Multiple Directionality options	•	•	•
Adaptation Manager	•	•	•
Fitting Formulas	VAC+, NAL-NL1/NAL-NL2, DSL 5.0	VAC+, NAL-NL1/NAL-NL2, DSL 5.0	VAC+, NAL-NL1/NAL-NL2, DSL 5.0
Connecting to the World			
Stereo streaming (2.4 GHz)	•	•	•
Oticon ON App & Oticon RemoteCare App	•	•	•
ConnectClip	•	•	•
EduMic	•	•	•
Remote Control 3.0	•	•	•
TV Adapter 3.0	•	•	•
Phone Adapter 2.0	•	•	•
Special Requirements			
Tinnitus SoundSupport™	•	•	•

Oticon miniRITE R: Oticon More™



oticon
life-changing technology

	1. Understanding speech 	2. Handling noise 	3. Unique fitting 	4. Handling environments 	5. Connectivity 
	<ul style="list-style-type: none"> • Speech clarity and richness • More adjustment options if you need it • Communication with facemasks 	<ul style="list-style-type: none"> • Contrast between speech and noise • Comfort in difficult situations 	<ul style="list-style-type: none"> • Customisation to your hearing fingerprint • Troubleshooting • Adjustments to changes in your hearing 	<ul style="list-style-type: none"> • Adjustments to fit the variability of your daily life environments 	<ul style="list-style-type: none"> • 2.4GHz • Made for iPhone & Android • ConnectClip • Options for TV, Telecoil & remote control
Oticon More 1	★★★★★★	★★★★★★	★★★★★★	★★★★★★	★★★★★★
Oticon More 2	★★★★★★	★★★★★★	★★★★★★	★★★★★★	★★★★★★
Oticon More 3	★★★★★★	★★★★★★	★★★★★★	★★★★★★	★★★★★★

In the full sound scene, Oticon More 1 can create greater contrast between meaningful sounds/speech and less meaningful sounds/noise than Oticon More 3.



Speech

Meaningful sounds

Less meaningful sounds/noise

